


ANAESTHETIC GAS SCAVENGING SYSTEM



Features

- GENTEC® Anaesthetic Gas Scavenging System is used for the safe removal of excess anaesthetic gas and anaesthetic vapors
- Color-coded for easy gas service identification according to ISO 32
- Complies with ISO 7396-2 and ISO 9170-2
- Flow: 50 l/min, 80 l/min (Pressure drop: 0.15~0.29 psi)
- Air driving flow : 10~15 l/min (Pressure: 58 psi)
- 100% tested and individually packaged

Dimensions (mm)



Ordering Information

Part Number	Model Number	Maximum Flow	Exhaust Diameter	Installation
A65216090	3816E-AGSS(1L)-HB16-B	50 L/min	I.D 16 mm Hose	Console Mount
A65216100	3816E-AGSS(1H)-HB16-B	80 L/min	I.D 16 mm Hose	Console Mount
A65226010	3826E-AGSS(1L)-HB22-W	50 L/min	I.D 22 mm Hose	Box Mount
A65226020	3826E-AGSS(1H)-HB22-W	80 L/min	I.D 22 mm Hose	Box Mount
A65226030	3826E-AGSS(1L)-HB16-F	50 L/min	I.D 16 mm Hose	Pannel Mount
A65226040	3826E-AGSS(1H)-HB16-F	80 L/min	I.D 16 mm Hose	Pannel Mount